

Bench mounted unit designed to demonstrate the use of ion-exchange resins for either continuous water softening or demineralisation.

The equipment is designed to show industrial operations such as 'breakthrough' and regeneration cycles.

UNIVERSAL SINGLE WIRE INTERFACE
DATA LOGGING AND CONTROL VIA PC
ACID PROOF BASE


W9MKII cation and anion resin in columns



armBUS USB port and sockets



armBUS Desktop Main Screen



Features / benefits

- ▶ Transparent test section for demonstration
- ▶ Control and data logging via PC
- ▶ Inline or outlet conductivity and pH sensors
- ▶ 5 x solenoid valves for automatic selection of each part of the process
- ▶ Operational manual with teaching exercises
- ▶ Universal single-wire interface

UK office - email: sales@armfield.co.uk tel: +44 (0) 1425 478781 (for ROW)
USA office - email: info@armfield.inc tel: +1 (609) 208-2800 (USA only)

Issue: 2
URL: <http://www.armfield.co.uk/w9>

Applications

ChE CE IP

We reserve the right to amend these specifications without prior notice. E&OE © 2019 Armfield Ltd. All Rights Reserved

Description

The unit consists of two vertical transparent columns containing the cation and anion resins, mounted on a support frame together with the necessary pumps, valves and sensors.

The arrangement of pumps and solenoid valves enables the flow configuration to be changed to demonstrate the various stages involved in the demineralisation or water-softening processes. Each stage of the process is selected by the operator using a PC and a mimic diagram for the appropriate stage displayed on the screen.

Conductivity, temperature and pH of the influent and effluent of the columns are measured using electronic sensors for data logging.

The apparatus is supplied with a typical commercial cation and anion resin, but other ion-exchange materials (not supplied) may be used so that their characteristics, exchange capacity etc. may be measured and compared.

Union couplings at both ends of the columns enable the column to be easily removed for maintenance and set-up without the need for tools. Regenerant, test and wash solutions are stored in separate containers to avoid cross-contamination, and effluents are fed to a waste-water tank for proper disposal.

Software

The ArmBUS software enables the operator to select the appropriate stage of the process and a mimic diagram with measured variables is displayed. The speed of the pump can be varied to meet the required flow rate.

Results are saved in a log, which can be viewed and manipulated with the ArmBUS results viewer. Results can be printed or exported in a spreadsheet format, which can be opened in a wide range of packages for further analysis.

Technical details

2 x vertical columns	16mm ID
Solenoid valves	5 x Two-way valves (suitable for 10% HCl, 5% NaOH and 10% NaCl)
Pumps	4 x 0-50mL (suitable for acid and alkali)
Pump	1 x 0-75mL (suitable for acid and alkali)
PH sensors	2
Conductivity/temperature sensors	2 x 0 – 2mS/cm
Conductivity probes	2 x 0 – 2mS/cm
Regenerant bottles	2 x 1L
Waste-water container	25L
Feed & wash water containers	2 x 10L

Overall dimensions

Length	0.90m
Width	0.45m
Height	0.80m
Packed and crated shipping specifications	
Volume	0.32m ³
Gross weight	35Kg

Requirements

Scale



- ▶ **Electrical supply:** supplied with a universal mains adaptor suitable for 100-240V / 1ph / 50-60Hz
- ▶ Software requires a computer running Windows XP or above with a USB port (computer not supplied by Armfield)
- ▶ Deionised water (not supplied)
- ▶ 5% sodium hydroxide (not supplied)
- ▶ 10% hydrochloric acid (not supplied)
- ▶ 10% sodium chloride (not supplied)

Water supply: Initial fill and drain



Learning objectives

- ▶ Water softening using cationic resin
- ▶ Exchange capacity of a cationic softening system
- ▶ Regeneration efficiency of a cationic softening system (regenerated using NaCl)
- ▶ Demineralisation using two-bed exchange (cationic resin and anionic resin used in series)
- ▶ Regeneration efficiency of a demineralisation system (regeneration cationic resin using HCl and anionic resin using NaOH)
- ▶ Exchange capacities of different resin materials (alternative materials not supplied)

Ordering specification

- ▶ Benchtop frame unit
- ▶ 2 x vertical transparent columns to house the resins
- ▶ 5 x peristaltic pumps for accurate process control and to prevent liquid contamination
- ▶ Two conductivity/temperature probes
- ▶ 2 x pH probes
- ▶ 5 x solenoid valves
- ▶ USB port for data logging using Armfield software
- ▶ Operational manual with teaching exercises
- ▶ Typical commercial cation and anion resins
- ▶ 2 x 1 L regenerant bottles
- ▶ 25 L waste-water container
- ▶ 2 x 10 L feed and wash water containers
- ▶ Powered by a universal power adaptor with worldwide approval

Ordering codes

- ▶ W9MklI

Armfield standard warranty applies with this product

Knowledge base

- > 28 years' expertise in research & development technology
- > 50 years' providing engaging engineering teaching equipment

Benefit from our experience, just call or email to discuss your laboratory needs, latest project or application.

An ISO 9001:2015 Company



armfield.co.uk

Aftercare

Installation
Commissioning
Training
Service and maintenance
Support: armfieldassist.com