

The apparatus consists of a rotating circular turntable attached onto the end of a vertical shaft. The shaft is held vertical within a sturdy wall bracket that must be attached to a rigid vertical surface, i.e. a wall or the SV100 Bench Mounted Frame.

At the lower tip of the shaft a special fixing has been created to accept a variety of end pivots. Mating seating's have been manufactured for the end pivots. The seating's positively locate at the bottom of the wall bracket so that they mate with their respective seating's but do not rotate under testing.

The seating's are made from different materials and with different internal cone angles.

Pure torsion is applied to the table through a pair of diametrically opposite loaded cords and pulleys. A set of calibrated weights and hangers are supplied.

To vary the thrust on the end pivot, the calibrated weights are added incrementally to the turntable surface. Weights are then added to the load hangers to calculate the necessary parameters.

#### Experimental content

- ▶ To investigate the relationship between frictional torque and axial thrust
- ▶ To determine the influence of the bearing cone angle
- ▶ To obtain the coefficient of friction for different bearings

#### Related laws

- ▶ Automotive
- ▶ Friction
- ▶ Machines
- ▶ Clutches
- ▶ Brakes
- ▶ Plate pressure
- ▶ Wear
- ▶ Torque
- ▶ Slip

#### Essential accessories / equipment

- ▶ SV100 Bench Mounted Frame

#### Features / benefits

- ▶ Wall mounted apparatus
- ▶ Can be mounted to the SV100 Bench Mounted Frame
- ▶ Pure torsion applied to the table
- ▶ Apparatus supplied with four sets of pivots
- ▶ The seating is in mild steel and the pivot is in brass
- ▶ Ball thrust bearing is included
- ▶ All load hangers and calibrated weights are supplied

#### Overall dimensions

Length	0.330m
Width	0.270m
Height	0.200m

#### Packed and crated shipping specifications

Volume	0.02m <sup>3</sup>
Gross weight	19kg

UK office - email: sales@armfield.co.uk tel: +44 (0) 1425 478781 (for ROW)

USA office - email: info@armfield.inc tel: +1 (609) 208-2800 (USA only)

A VERTICALLY MOUNTED EXPERIMENT TO INVESTIGATE THE RELATIONSHIP BETWEEN FRICTIONAL TORQUE AND AXIAL THRUST



#### Requirements

SV  
100

#### Scale



- ▶ Rigid vertical support

#### Technical specification

- ▶ Four pivot sets at angles of 60°, 90°, 120°, 180° (flat)
- ▶ Turntable mass: Approximately 1kg
- ▶ Turntable effective diameter: Ø200mm
- ▶ Thrust bearing seating supplied

#### Ordering specification

- ▶ SV902 Pivot Friction Apparatus
- ▶ 4 x Pivot set
- ▶ 2 x Load hanger
- ▶ 1 x Thrust bearing set
- ▶ 1 x 10m Spare cord
- ▶ 10 x 0.1N; 10 x 0.2N; 2 x 1N; 8 x 2N; 2 x 50N; 1 x 100N
- ▶ Instruction manual

#### Ordering codes

- ▶ SV902 Pivot Friction Apparatus
- ▶ SV100 Bench Mounted Frame

Issue: 1

URL: <http://www.armfield.co.uk/structures>

Applications

ME CE IP

We reserve the right to amend these specifications without prior notice. E&OE © 2022 Armfield Ltd. All Rights Reserved